



**GELLER
MICROANALYTICAL
LABORATORY, Inc.**

426e BOSTON ST., TOPSFIELD, MA 01983-1216
TEL 978 887-7000 FAX 978-887-6671

Jg@gellermicro.com <http://www.gellermicro.com>

Certified to ISO-9001 and 17025



VUS-J-XL

Extra Large Sized Vacu-Storr™

Measuring 6" diameter X 12" high this Vacu-Storr holds vacuum for a long long time!

The photo at the left shows the Vacu-Storr with protective netting (provided standard) and a polymer coating over the glass. The photo below shows a Vacu-Storr that was dropped 6 ft. onto pavement. The vessel fractured but almost all the glass fragments are captured by the netting covered polymer coating.

